

BLOCKING MATERIALS

BOND365®



Bond365® Part A

Bond365® Part B



Uni-bloc Wood Block



Toe Cup



Toe Cup Holder



Mixing Sticks

Helpful Hint:

Cold out? Keep Part B in your pocket and Part A in a warming bag prior to application. **Bond365®** performs best between 65 - 75° F.

APPLICATION INSTRUCTIONS

STEP ONE



Dispense entire contents of Part B into toe cup.

STEP TWO



Dispense entire contents of Part A into toe cup.

STEP THREE



Mix the two parts vigorously with mixing stick for 25 seconds. **Bond365®** will turn opaque in color.

STEP FOUR



When **Bond365®** becomes clear, insert block into bottom of toe cup.

STEP FIVE



Apply thin layer to both block and claw. Spread remaining glue on block to ensure 3/16—1/4" (5-7 mm) thickness throughout.

STEP SIX



Position block on prepared claw. Hold until **Bond365®** begins to harden. Reposition the block as needed. Allow to cure for four (4) minutes before weight bearing.